HeyTec Flow-V

(hochviskoses Flowable Composite)

- (Highly Viscous Flowable Composite)
- (composite fluide hautement visqueux)
- (composito flowable altamente viscoso)
- (composite fluido, altamente viscosa)
- (Высоковязкий текучий композит)

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- Fissuren sealing
- Extensive fissure sealing on molars and premolars
- Prevention of class V cervical caries, root erosions, v-shaped lesions
- Minimally invasive fillings in Black's classes I, II and III
- Corrections of enamel defects
- Blocking out undercuts
- Small corrections of shape and color of the enamel

Procedure:

Pre-treatment:
Before starting the treatment, clean the hard tooth substance with a fluoride-free polishing paste. Select the correct color with help of the Vita® shade guide while the tooth is still moist.

1. Cavity preparation

Prepare the cavity according to the general rules of the adhesive technique. Use a procedure that is gentle on the hard tooth substance. Bevel the incisal edges on anterior teeth. Do not bevel the edges on posterior teeth. Avoid irregular edges. Next, clean the cavity with water spray, remove all residue and dry. It is necessary to dry the site completely. We recommend the use of a kofferdam.

2. Pulp protection/liner

The use of a liner can be foregone, if an enamel dentin adhesive is used. In case of very deep cavities that come close to the pulp, cover the according surfaces with a calcium hydroxide compound.

3. Design of the approximal contacts

For cavities in approximal areas, insert and fix a transparent matrix.

4. Adhesive system

Etch and bond according to the manufacturer's instructions.

5a. Application of HeyTec Flow-V syringes

Fill HeyTec Flow-V in thin layers (max. 2 mm) directly into the cavity. Use the dispensing curved application tips included with the material. For hygiene reasons, do not use the application tips more than one time. Avoid introducing air bubbles into the material during application with the tip. Ensure that the prepared hard tooth substance is well-coated.

5b. Application of HeyTec Flow-V compules

Insert the composite into the dispenser. Remove the sealing cap. Place the composite into the cavity and move it towards the cavity. Press the material into the cavity. Apply slow and steady pressure to the compule. Do not use excessive force! (maximum layer thickness 2 mm). To remove the compule from the dispenser after use, retract the plunger. Next, remove the compule. Please note: For hygiene reasons, compules are intended for single use only.

6. Polymerization

Polymerizationsschicht layer with a conventional polymerization lamp for 40 sec. Hold the light guide as close as possible to the surface of the filling. A dispersion layer will form on the surface of the material during polymerization. This dispersion layer forms the chemical bond between the layers and must not be touched or removed.

7. Finishing

HeyTec Flow-V is ready to be finished and polished after curing. Suitable tools are finishing diamonds, flexible wheels, silicone polishers and polish-ing brushes.

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- In case of extensive restorations, the surgical light should be moved away from the working area temporarily to avoid premature curing of the composite. Alternatively, the material can be covered with a light-tight foil.
- For hygiene reasons, do not use the curved application tips included with the product more than once.
- Use the light guide with the emission range of 350–500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions.

Light intensity for polymerization ≥ 650 mW/cm²

Wavelength for polymerization 350–500 nm

Polymerization time 40 sec.

Hazard statement / precautionary statement:

Contains tetramethylene dimethacrylate

Warning: May cause an allergic skin reaction. Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

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Instructions for use

HeyTec Flow-V is a light-curing, flowable, highly radiopaque (210% Al) composite with a high viscosity. The guidelines and requirements of DIN EN ISO 4049 apply. HeyTec Flow-V is available in syringes and compules. The compules are for single use only. Please do not reuse them, as this makes it impossible to rule out contamination and germ infection.

In case of exposed dentin in a cavity, ensure a proper protection of the pulp to avoid possible pulpal reactions (e.g. g. line the cavity with a calcium hydroxide compound).

Composition:

Glass powder, diurethane dimethacrylate, tetramethylene dimethacrylate, silicon dioxide

Total filler:

77 % by weight (57 % by volume) of inorganic filler (0.005–40 µm)

Indications:

- Fissure sealing
- Extensive fissure sealing on molars and premolars
- Prevention of class V cervical caries, root erosions, v-shaped lesions
- Minimally invasive fillings in Black's classes I, II and III
- Corrections of enamel defects
- Blocking out undercuts
- Small corrections of shape and color of the enamel

Procedure:

Pre-treatment:
Before starting the treatment, clean the hard tooth substance with a fluoride-free polishing paste. Select the correct color with help of the Vita® shade guide while the tooth is still moist.

1. Cavity preparation

Prepare the cavity according to the general rules of the adhesive technique. Use a procedure that is gentle on the hard tooth substance. Bevel the incisal edges on anterior teeth. Do not bevel the edges on posterior teeth. Avoid irregular edges. Next, clean the cavity with water spray, remove all residue and dry. It is necessary to dry the site completely. We recommend the use of a kofferdam.

2. Pulp protection/liner

The use of a liner can be foregone, if an enamel dentin adhesive is used. In case of very deep cavities that come close to the pulp, cover the according surfaces with a calcium hydroxide compound.

3. Design of the approximal contacts

For cavities in approximal areas, insert and fix a transparent matrix.

4. Adhesive system

Etch and bond according to the manufacturer's instructions.

5a. Application of HeyTec Flow-V syringes

Fill HeyTec Flow-V in thin layers (max. 2 mm) directly into the cavity. Use the dispensing curved application tips included with the material. For hygiene reasons, do not use the application tips more than one time. Avoid introducing air bubbles into the material during application with the tip. Ensure that the prepared hard tooth substance is well-coated.

5b. Application of HeyTec Flow-V compules

Insert the composite into the dispenser. Remove the sealing cap. Place the composite into the cavity and move it towards the cavity. Press the material into the cavity. Apply slow and steady pressure to the compule. Do not use excessive force! (maximum layer thickness 2 mm). To remove the compule from the dispenser after use, retract the plunger. Next, remove the compule. Please note: For hygiene reasons, compules are intended for single use only.

6. Polymerization

Polymerizationsschicht layer with a conventional polymerization lamp for 40 sec. Hold the light guide as close as possible to the surface of the filling. A dispersion layer will form on the surface of the material during polymerization. This dispersion layer forms the chemical bond between the layers and must not be touched or removed.

7. Finishing

HeyTec Flow-V is ready to be finished and polished after curing. Suitable tools are finishing diamonds, flexible wheels, silicone polishers and polish-ing brushes.

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Instructions for use

HeyTec Flow-V is a light-curing, flowable, highly radiopaque (210% Al) composite with a high viscosity. The guidelines and requirements of DIN EN ISO 4049 apply. HeyTec Flow-V is available in syringes and compules. The compules are for single use only. Please do not reuse them, as this makes it impossible to rule out contamination and germ infection.

In case of exposed dentin in a cavity, ensure a proper protection of the pulp to avoid possible pulpal reactions (e.g. g. line the cavity with a calcium hydroxide compound).

Composition:

Phenolic substances (e.g. eugenol) inhibit polymerization.

Therefore, these materials (e.g. zinc oxide eugenol cements) must not be used as liners.

Disinfection / Protection from cross-contamination:

Place the functioning syringe with attached application tip into a suitably shaped barrier sheath; pierce end of sheath with cannula, enclosing the cannula for disinfection. Using a barrier sheath facilitates cleaning and disinfection of the syringe between patients.

After use of sheathed syringe, remove delivery tip and sheath by grasping on the hub of the delivery tip through the sheath, twist and remove tip along with sheath. Discard tip and sheath in appropriate waste stream. Replace syringe storage cap.

Disinfect – After removal of the delivery tip and sheath, disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Health-Care Settings – 2003 (Vol. 52, No. RR-17), Centers for Disease Control and Prevention (USA).

Storage:

Store at temperatures between 10 °C and 25 °C (50 °F to 77 °F). Screw the cap back onto the syringe tightly after each use. Let the material reach room temperature before use. Withdraw the plunger slightly after use to keep the outlet from becoming plugged. Do not use after the expiration date (see label on the syringe). Only for use in dentistry. Keep out of children's reach. This product has been developed for the specific use illustrated above. Only process as described in these instructions. The manufacturer will not be held liable for any damages that result from improper use or improper processing.

5a. Application of HeyTec Flow-V syringes

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*Vita es una marca registrada de Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.